



LogTag Recalibration

Florida Department of Health

Immunization Section

Field Operations

Lori Wright, Area 8 Field Office Manager

June 22, 2016





LogTag Recalibration

Step 1- Submit





Time to Recalibrate?





Time to Recalibrate?

CALIBRATION CERTIFICATE LOG

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS
Freezer: 291010-REFRIGERATOR

Selection criteria

Expiration Date From: 06/11/2011 x Thru: 06/08/2021

Get matching certificates

Current For Unit	Expiration Date	Accredited Lab	Uncertainty	Pass/Fail	----- Non-accredited Lab Only -----	ISO 17025	
<input checked="" type="checkbox"/>	11/7/2016	PJLA	Yes		----- n/a -----		Delete
<input type="checkbox"/>			Yes		----- n/a -----		Delete

Return to Refrigerators/Freezers List



Certificate of Calibration

Report Number: 14-1107-01
Calibration Date: Friday, November 7, 2014
Calibration Due Date: November 7, 2016

Customer: Florida Department of Health
Customer PO #: AB3948
Sales Order #: 60802

Manufacturer: LogTag
Model: TRED30-7R
Serial Number: 1060019051



Time to Recalibrate?

Certificate of calibration for the thermometer attached to storage unit PIN:700140-SPIDER2.1-FREEZERA2 will expire on 04/06/2016. You must recalibrate your thermometer BEFORE or BY the calibration expiry date to continue ordering and receiving VFC vaccines. Please contact your Area Immunization Consultant if you need additional details on recalibrating your thermometer.

Please be advised that another Calibration record is linked to this unit which expires on 07/05/2016.

Thanks
Florida VFC Program

I



Time to Recalibrate?

Page 1 of 1



879 Maple St. Hopkinton, New Hampshire 03229



Accreditation # 76289

Certificate of Calibration

Report Number: 14-1107-01

Calibration Date: Friday, November 7, 2014

Calibration Due Date: November 7, 2016

Customer: Florida Department of Health

Customer PO #: AB3948

Sales Order #: 60802

Manufacturer: LogTag

Model: TRED30-7R

Serial Number: 1060019051

Example: Calibration expires 11/7/16. Submit for recalibration no sooner than 9/7/16 (60 days).

How many LogTags Can I Return at Once?



One Thermometer with 2 probes= 2 thermometers.
You can return 2 LogTags at once.

Or, two thermometers with a
single probe = 2 thermometers.
You can return 2 LogTags at
once.



How many LogTags Can I Return at Once?



One Thermometer with one probe. If this is the only backup thermometer, return one LogTag at a time.



Begin the Recalibration Process



Send

To... VFCRecalibration@flhealth.gov

Cc...

Bcc...

Subject: 290010 Recalibration

Begin the Recalibration Process



Send

To... VFCRecalibration@flhealth.gov

Cc...

Bcc...

Subject: 290010 Recalibration

Refrigerator Log Tag 290010-REFRIGERATOR s/n 1360001212 and Freezer Log Tag s/n 1360002212 need recalibration.

Begin the Recalibration Process



From: zzzz Feedback, HSD_FloridaVFCRecalibration
To: Osgood, Norma
Cc:
Subject: RE: PIN 700069 Recalibration

Sent: Tue 6/7/2016 1:24 PM

Thank you for contacting the Immunization Section. You are receiving this email in response to your request to recalibrate one or more of your LogTag data loggers.

To begin this process:

- Click on the link below to be redirected to the MicroDAQ website.
- Complete the online form. If you have two currently certified backup thermometers, both LogTags for your PIN can be included on the same form, just be sure to click Yes on "calibrate additional logger", and type the second serial number in the comments section. (If there is only one currently certified backup thermometer, you will have to send one now, and one when the first one is returned.)
- Click the "submit" button.
- MicroDAQ will confirm the receipt of your form and email further instructions for mailing.

<https://www.microdaq.com/fl-calibration>

As a reminder, only LogTags supplied by the VFC Program are eligible for this service. You, as the provider, are responsible only for the initial shipping charges. All other costs will be paid by VFC.

Please do not reply to this email. If you have any further questions, please contact your local Field Staff Representative.

VFC Recalibration Team

Begin the Recalibration Process



From: zzzz Feedback, HSD_FloridaVFCRecalibration
To: Osgood, Norma
Cc:
Subject: RE: PIN 700069 Recalibration


Sent: Tue 6/7/2016 1:24 PM

Thank you for contacting the Immunization Section. You are receiving this email in response to your request to recalibrate one or more of your LogTag data loggers.

To begin this process:

- Click on the link below to be redirected to the MicroDAQ website.
- Complete the online form. If you have two currently certified backup thermometers, both LogTags for your PIN can be included on the same form, just be sure to click Yes on "calibrate additional logger", and type the second serial number in the comments section. (If there is only one currently certified backup thermometer, you will have to send one now, and one when the first one is returned.)
- Click the "submit" button.
- MicroDAQ will confirm the receipt of your form and email further instructions for mailing.

<https://www.microdaq.com/fl-calibration>



As a reminder, only LogTags supplied by the VFC Program are eligible for this service. You, as the provider, are responsible only for the initial shipping charges. All other costs will be paid by VFC.

Please do not reply to this email. If you have any further questions, please contact your local Field Staff Representative.

VFC Recalibration Team

Begin the Recalibration Process



General Information

i Contact Name*

Type Your Full Name

i Practice Name*

Enter Your Practice Name

i Address*

Enter Your Address

Address

i City*

Enter Your City

i State

FL

i Zip Code*

Enter Your Zip Code

i Telephone*

Enter Your Telephone Number

i E-Mail*

Enter Your Email Address: example@mail.

Calibration Information

Calibration Type

Fridge



Serial Number

Enter the Serial Number Listed on Device

Calibrate Additional Data Logger

☒ No

☐ Yes

Comments

i Comments

Enter Any Additional Information

Begin the Recalibration Process



Calibration Information

Calibration Type

Fridge

Serial Number

Enter the Serial Number Listed on Device

Calibrate Additional Data Logger

☒ No

☐ Yes

Comments

Comments

Enter Any Additional Information

Begin the Recalibration Process



Submit

Begin the Recalibration Process



MicroDAQ and MDQ Calibration Lab has received your FL DOH Calibration Order for Your TRED30 Vaccine Data Logger

Print this email and include a copy in the Box with the Data Logger and Probe/Vial.

Typical Lead Time for Recalibration: 7 to 10 Business Days after Receipt of Data Loggers

Instructions:


Ship the Logger Back to Us Using a Trackable Service (UPS, FedEx, etc.)

You Only Need to Send Us the Logger and the Probe/Vial. We do NOT Need the Docking Station or Any Other Accessories


Ship the package, using a trackable service (UPS or FedEx) to:

MDQ Calibration Lab
Ref: FL DOH
879 Maple Street
Contoocook, NH 03229

Begin the Recalibration Process



Begin the Recalibration Process



A new LogTag has been detected in the interface. Please wait while its readings are downloaded.

There are 1 remaining LogTag(s) to download.

Cancel

Begin the Recalibration Process



Temperature Log for Refrigerator (Fahrenheit)

Day 1 – 15 Month/Year: 9/16
 VFC PIN #: 291010

- Completing this temperature log:** Check the temperature in the refrigerator compartment of your vaccine storage unit at least twice each working day. Place an "X" in the box that corresponds with the temperature, the time of the temperature reading, and your initials. Save each month's completed form for 3 years, unless state or local jurisdictions require a longer time period.
- Temperatures warmer than 46°F or colder than 35°F are outside the acceptable temperature range. You must take action!**

Take Action!
 If temperature is **too warm** (above 46°F) or **too cold** (below 35°F):

- Store the vaccine** under proper conditions as quickly as possible.
- Temporarily mark exposed vaccine "do not use"** until you have verified whether or not the vaccine may be used.
- Call the VFC Program** for further assistance at 1-800-483-2543.
- Document the action taken** on the reverse side of this log.

VFC Customer Service Representative

Immunization Section Area Field Staff Office

Staff Initials	[Handwritten: Z, Z, Z, Z, Z]					[Handwritten: RK, LF, RR, RR, RR]									
Day of Month	[Handwritten: 6, 7, 8, 9]														
Room Temp	[Handwritten: 45, 45, 45, 45]														
Exact Time	[Handwritten: 8:00 am, 4:30 pm, 8:00 am, 4:30 pm, 7:30 pm]														
Min/Max Temp (since previous reading)	[Handwritten: 45, 45, 45, 45]														

*Write any unacceptable temps (above 46°F or below 35°F) on the lines below. Then take action!

Danger! Temperatures above 46°F are too warm! Write any unacceptable temperature on the lines above* and call your local VFC Program immediately!

Acceptable Temperatures	46°															
	45°															
	44°															
	43°															
	42°															
	41°															
	40°															
	39°															
	38°															
	37°															
	36°															
	35°															

Danger! Temperatures below 35°F are too cold! Write any unacceptable temperature on the lines above* and call your local VFC Program immediately!

See back for "Vaccine Storage Troubleshooting Record"

Vaccines for Children
Protecting All of Florida's Children

Florida Vaccines for Children Program
 4052 Bald Cypress Way, Bin A-11
 Tallahassee, FL 32399

Phone: 1-800-483-2543, Fax: 850-245-4734
 Email: FloridaVFC@FLHealth.gov
www.immunizeflorida.org/vfc

Florida HEALTH RF1-15

Package the LogTags



Instructions:

Ship the Logger Back to Us Using a Trackable Service (UPS, FedEx, etc.)

You Only Need to Send Us the Logger and the Probe/Vial. We do NOT Need the Docking Station or Any Other Accessories

Ship the package, using a trackable service (UPS or FedEx) to:

MDQ Calibration Lab
Ref: FL DOH
879 Maple Street
Contoocook, NH 03229

Shipping Instructions

1. Use a Sturdy Box - NO Packing Envelopes 2. Disconnect Probe from Sensor Before Shipping 3. Insert Glycol Bottle and Probe Cable into a Zip-Lock Bag 4. Use a Rubber Band to Secure the Logger to the Outside of the Sealed Bag 5. Wrap the Logger and Bag with Packing Materials (Newspaper, Packing Peanuts, etc.) to Protect them in the Box

Package the LogTags



Package the LogTags



Package the LogTags



Package the LogTags



Package the LogTags



Instructions:

Ship the Logger Back to Us Using a Trackable Service (UPS, FedEx, etc.)

You Only Need to Send Us the Logger and the Probe/Vial. We do NOT Need the Docking Station or Any Other Accessories

Ship the package, using a trackable service (UPS or FedEx) to:

MDQ Calibration Lab
Ref: FL DOH
879 Maple Street
Contoocook, NH 03229

Shipping Instructions

1. Use a Sturdy Box - NO Packing Envelopes 2. Disconnect Probe from Sensor Before Shipping 3. Insert Glycol Bottle and Probe Cable into a Zip-Lock Bag> 4. Use a Rubber Band to Secure the Logger to the Outside of the Sealed Bag 5. Wrap the Logger and Bag with Packing Materials (Newspaper, Packing Peanuts, etc.) to Protect them in the Box


Package the LogTags



*Turnaround time for MicroDAQ to recalibrate your LogTag(s) and return them to you is 7-10 business days.

Finishing



 Send

To... VFCRecalibration@flhealth.gov

Cc...

Bcc...

Subject: 290010 Recalibration

291010-REFRIGERATOR has been shipped to MicroDAQ for recalibration.

Finishing



CALIBRATION CERTIFICATE LOG

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS

Freezer: 291010-REFRIGERATOR

Selection criteria

Expiration Date From: 06/11/2011 × Thru: 06/08/2021

Get matching certificates

Current For Unit	Expiration Date	Accredited Lab	Uncertainty	Pass/Fail	Non-accredited Lab Only	ISO 17025	
<input checked="" type="checkbox"/>	11/7/2016	PJLA	Yes		----- n/a -----		Delete
<input type="checkbox"/>	02/20/2017	PJLA	Yes		----- n/a -----		Delete

Return to Refrigerators/Freezers List

Finishing




REFRIGERATORS / FREEZERS LIST

[Show Help Text](#)

VFC Provider: * XYZ PEDIATRICS

☒ Pin ☐ VFC EID

PIN/EID	Identifier	R/F	Ended Date	<input type="checkbox"/> Show out-of-service units	
P: 291010	291010-REFRIGERATOR	R			Upload Temperatures File Add Temp Reading Add Missing Temp Reading Temperature Log Excursion Incident List Add Calibration Calibration Log !
				<input type="button" value="Show Detail"/>	
P: 291010	291010-FREEZER	F			Upload Temperatures File Add Temp Reading Add Missing Temp Reading Temperature Log Excursion Incident List Add Calibration Calibration Log !
				<input type="button" value="Show Detail"/>	



Finishing



CALIBRATION CERTIFICATE LOG

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS

Freezer: 291010-REFRIGERATOR

Selection criteria

Expiration Date From: Thru: [Get matching certificates](#)

Current For Unit	Expiration Date	Accredited Lab	Uncertainty	Pass/Fail	Non-accredited Lab Only	ISO 17025	
<input checked="" type="checkbox"/>	11/7/2016	PJLA	Yes		----- n/a -----		Delete
<input type="checkbox"/>	02/20/2017	PJLA	Yes		----- n/a -----		Delete

[Return to Refrigerators/Freezers List](#)

Finishing



Fix exp date to 7/11/16

UPDATE CALIBRATION RECORD

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS
Refrigerator: 291010-REFRIGERATOR

Help with completing this form: [Calibration Report Examples](#)

Thermometer

Manufacturer/Name: LOGTAG
Model Nbr: * TRED30-7R
Serial Nbr: * 1360008146

Expiration Date: * 01/07/2018 x

Is Current for the unit: ☒

☒ Certificate is from an accredited laboratory

Laboratory: * PJLA v

☒ Documented Uncertainty.

Measurement results indicate unit passed test and the documented uncertainty is within suitable limits (recommended uncertainty = +/- 1F (0.5C)).

[Change](#)

[Return to Calibration Log](#)

Finishing



CALIBRATION CERTIFICATE LOG

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS

Freezer: 291010-REFRIGERATOR

Selection criteria

Expiration Date From: 06/09/2011 x

Thru: 06/06/2021

Get matching certificates

Current For Unit	Expiration Date	Accredited Lab	Uncertainty	Pass/Fail	----- Non-accredited Lab Only ----- ISO 17025
<input type="checkbox"/>	11/7/2016	PJLA	Yes		----- n/a -----
<input checked="" type="checkbox"/>	02/20/2017	PJLA	Yes		----- n/a -----

Return to Refrigerators/Freezers List

What If My Backup Thermometer Isn't Listed?



CALIBRATION CERTIFICATE LOG

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS
Freezer: 291010-REFRIGERATOR

[Selection criteria](#)

Expiration Date From: 06/11/2011

Thru: 06/08/2021

Current For Unit	Expiration Date	Accredited Lab	Uncertainty	Pass/Fail	----- Non-accredited Lab Only ----- ISO 17025	
<input checked="" type="checkbox"/>	11/7/2016	PJLA	Yes		----- n/a -----	<input type="button" value="Delete"/>

What If My Backup Thermometer Isn't Listed?



REFRIGERATORS / FREEZERS LIST

[Show Help Text](#)

VFC Provider: * XYZ PEDIATRICS

☒ Pin ☐ VFC EID

[Get list for Pin](#)

PIN/EID	Identifier	R/F	Ended Date	<input type="checkbox"/> Show out-of-service units	
P: 291010	291010-REFRIGERATOR	R			Upload Temperatures File Add Temp Reading Add Missing Temp Reading Temperature Log Excursion Incident List Add Calibration Calibration Log !
				Show Detail	
P: 291010	291010-FREEZER	F			Upload Temperatures File Add Temp Reading Add Missing Temp Reading Temperature Log Excursion Incident List Add Calibration Calibration Log !
				Show Detail	



Review

- * For this step in the recalibration process:
 - * Within 30-60 days of your LogTag's Calibration Due Date, email vfc recalibration@flhealth.gov and wait for our reply
 - * The MicroDAQ website will be provided in our return email. Click the link in this email to go to the MicroDAQ website
 - * Complete the online form
 - * Print the confirmation email from MicroDAQ
 - * Download the data from the LogTag and upload to Florida SHOTS
 - * Place the backup thermometer into the storage unit (If you have not added your backup thermometer to this storage unit's calibration log, please do so now. Do not delete your Log Tag)
 - * Ship your LogTags to MicroDAQ, packaged as described



Review

- * For this step in the recalibration process (cont.):
 - * Email vfc recalibration@flhealth.gov to let us know which LogTag(s) are being recalibrated
 - * Manually document temperatures twice daily
 - * Make your backup thermometer “current for unit” in Florida SHOTS
 - * Enter temperatures at least weekly until your LogTag(s) are returned



Step 2- LogTags Are Returned





Step 2- LogTags Are Returned





Step 2- LogTags Are Returned





Step 2- LogTags Are Returned

REFRIGERATORS / FREEZERS LIST

[Show Help Text](#)

VFC Provider: * XYZ PEDIATRICS

☒ Pin ☐ VFC EID

[Get list for Pin](#)

PIN/EID	Identifier	R/F	Ended Date	<input type="checkbox"/> Show out-of-service units	
P: 291010	291010-REFRIGERATOR	R			Upload Temperatures File Add Temp Reading Add Missing Temp Reading Temperature Log Excursion Incident List Add Calibration Calibration Log !
				Show Detail	
P: 291010	291010-FREEZER	F			Upload Temperatures File Add Temp Reading Add Missing Temp Reading Temperature Log Excursion Incident List Add Calibration Calibration Log !
				Show Detail	



Step 2- LogTags Are Returned

REFRIGERATORS / FREEZERS LIST

[Show Help Text](#)

VFC Provider: * XYZ PEDIATRICS

☒ Pin ☐ VFC EID

PIN/EID	Identifier	R/F	Ended Date	<input type="checkbox"/> Show out-of-service units	
P: 291010	291010-REFRIGERATOR	R			Upload Temperatures File Add Temp Reading Add Missing Temp Reading Temperature Log Excursion Incident List Add Calibration Calibration Log !
				<input type="button" value="Show Detail"/>	
P: 291010	291010-FREEZER	F			Upload Temperatures File Add Temp Reading Add Missing Temp Reading Temperature Log Excursion Incident List Add Calibration Calibration Log !
				<input type="button" value="Show Detail"/>	



Step 2- LogTags Are Returned

CALIBRATION CERTIFICATE LOG

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS

Freezer: 291010-REFRIGERATOR

Selection criteria


Expiration Date From: 06/11/2011 x

Thru: 06/08/2021

Get matching certificates

Current For Unit	Expiration Date	Accredited Lab	Uncertainty	----- Non-accredited Lab Only -----	
				Pass/Fail	ISO 17025
<input checked="" type="checkbox"/>	11/7/2016	PJLA	Yes	----- n/a -----	Delete
<input type="checkbox"/>	02/20/2017	PJLA	Yes	----- n/a -----	Delete

Return to Refrigerators/Freezers List





Step 2- LogTags Are Returned

Provider Pin: 291010 Site XYZ PEDIATRICS

Freezer: 291010-REFRIGERATOR

Help with completing this form: [Calibration Report Examples](#)

Thermometer

Manufacturer/Name: LOGTAG
Model Nbr: * TRED30-7R
Serial Nbr: * 1360004338

Expiration Date: * 11/07/2016

Is Current for the unit: ☐

☒ Certificate is from an accredited laboratory

Laboratory: * PJLA

☒ Documented Uncertainty.

Measurement results indicate unit passed test and the documented uncertainty is within suitable limits (recommended uncertainty = +/- 1F (0.5C).



Step 2- LogTags Are Returned

UPDATE CALIBRATION RECORD

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS

Freezer: 291010-REFRIGERATOR

Help with completing this form: [Calibration Report Examples](#)

Thermometer

Manufacturer/Name: LOG TAG

Model Nbr: * TRED30-7R

Serial Nbr: * 1060016692

Expiration Date: * 11/07/2016

Is Current for the unit: ☐

☒ Certificate is from an accredited laboratory

Laboratory: * PJLA

☒ Documented Uncertainty.

Measurement results indicate unit passed test and the documented uncertainty is within suitable limits (recommended uncertainty = +/- 1F (0.5C).

Change

Return to Calibration Log

Cancel



Step 2- LogTags Are Returned

Page 1 of 1



879 Maple St. Hopkinton, New Hampshire 03229



Certificate of Calibration

Report Number: 14-1107-01

Generation Date: Friday, November 7, 2014

Calibration Due Date: November 7, 2018

Customer: Florida Department of Health

Customer PO #: AB3948

Sales Order #: 60802

Manufacturer: LogTag

Model: TRED30-7R

Serial Number: 1060019051



Step 2- LogTags Are Returned

Provider Pin: 291010 Site XYZ PEDIATRICS

Freezer: 291010-REFRIGERATOR

Help with completing this form: [Calibration Report Examples](#)

Thermometer

Manufacturer/Name:	LOGTAG
Model Nbr: *	TRED30-7R
Serial Nbr: *	1360004338

Expiration Date: * 11/07/2018

Is Current for the unit: ☐

☒ Certificate is from an accredited laboratory

Laboratory: * PJLA

☒ Documented Uncertainty.

Measurement results indicate unit passed test and the documented uncertainty is within suitable limits (recommended uncertainty = +/- 1F (0.5C).



Step 2- LogTags Are Returned

Provider Pin: 291010 Site XYZ PEDIATRICS

Freezer: 291010-REFRIGERATOR

Help with completing this form: [Calibration Report Examples](#)

Thermometer

Manufacturer/Name:	LOGTAG
Model Nbr: *	TRED30-7R
Serial Nbr: *	1360004338

Expiration Date: * 11/07/2018

Is Current for the unit: ☐

☒ Certificate is from an accredited laboratory

Laboratory: * PJLA

☒ Documented Uncertainty.

Measurement results indicate unit passed test and the documented uncertainty is within suitable limits (recommended uncertainty = +/- 1F (0.5C).

[Change](#)

[Return to Calibration Log](#)

[Cancel](#)



Step 2- LogTags Are Returned

CALIBRATION CERTIFICATE LOG

[Show Help Text](#)

Provider Pin: 291010 Site XYZ PEDIATRICS
Freezer: 291010-REFRIGERATOR

[Selection criteria](#)


Expiration Date From: 06/12/2011 x Thru: 06/09/2021 [Get matching certificates](#)

Current For Unit	Expiration Date	Accredited Lab	Uncertainty	Pass/Fail	Non-accredited Lab Only	ISO 17025	
<input type="checkbox"/>	11/7/2018	PJLA	Yes		----- n/a -----		Delete
<input checked="" type="checkbox"/>	2/20/2017	PJLA	Yes		----- n/a -----		Delete

[Return to Refrigerators/Freezers List](#)




Step 2- LogTags Are Returned

 Send	To...	VFCRecalibration@flhealth.gov
	Cc...	
	Bcc...	
Subject:		PIN 291010 Recalibrated

Refrigerator Log Tag 291010-REFRIGERATOR, serial number 1360001212 was received from MicroDAQ.

|





Step 2- LogTags Are Returned

LogTag

User ID: 291010-REFRIGERATOR

Push button start LogTag battery: OK

☒ Enable pre-start logging

☐ Record readings continuously, overwrite oldest when memory full

☒ Record readings so that:

Readings recorded will span at least 81 days

Number of readings to record 7,777 maximum is 7777

Record a reading every 15 Minutes

Begin recording after a delay of 5 Minutes

☒ Trigger alarm when readings above/equal 5 °F

☐ After 2 Accumulative violation readings (None)

☒ Trigger alarm when readings below/equal -40 °F

☐ After 2 Accumulative violation readings (None)

☐ Configure requires a password

☐ Download requires a password


Advanced Options

< Back Next > Close Help

REFRIGERATORS / FREEZERS LIST	
VFC Provider: * 1291010 XYZ PEDIATRICS	
<input checked="" type="radio"/> Pin <input type="radio"/> VFC EID	Get list for Pin
PIN/EID	Identifier
P: 291010	291010-REFRIGERATOR
P: 291010	291010-FREEZER




Step 2- LogTags Are Returned






Step 2- LogTags Are Returned






Step 2- LogTags Are Returned






Step 2- LogTags Are Returned






Step 2- LogTags Are Returned






Step 2- LogTags Are Returned






Step 2- LogTags Are Returned






Step 2- LogTags Are Returned






Step 2- LogTags Are Returned







Step 2- LogTags Are Returned






Step 2- LogTags Are Returned





Step 2- LogTags Are Returned





Step 2- LogTags Are Returned





Review For Part 2

- * For this step in the recalibration process:
 - * Check the LogTag and its probe for any sign of damage
 - * Place the probe into the appropriate storage unit
 - * Update the LogTag's new expiration date
 - * Make the LogTag "current for unit"
 - * Send an email to VFCRecalibration@flhealth.gov to let us know that this LogTag has been received. Please include your PIN and the unit ID.
 - * Reconfigure the logger and check the settings
 - * Wait 2-3 hours before starting the LogTag to allow the probe to acclimate to the storage unit



Questions or Assistance

- * If you have any questions about this process, please email vfc recalibration@flhealth.gov
- * PLEASE INCLUDE YOUR VFC PIN AND THE UNIT ID FOR THE LOGTAG IN QUESTION ON ALL CORRESPONDENCE WITH THE RECALIBRATION TEAM!